


GE Medical Systems

MODALITY INSTALLATION PLANNING 10/18/00, REVISION 8

SECTION 9 PAGE 9-10

SENOGRAPHE DMR



ROOM DIMENSIONS	LENGTH & WIDTH	CEILING HEIGHT
MINIMUM:	10'-0" x 12'-0" [3.5m x 3.7m]	7'-8" [2.3m]



	ROOM PLAN LOCATED ON PAGE 9-10						
ITEM	ITEM QUANTITY						
	†	ITEM DESCRIPTION	WEIGHT	HEAT OUTPUT			
1 2 3 4 50	1 1 1 1	EDM TOWER FOR SENOGRAPHE DMR SENOGRAPHE DMR GENERATOR CABINET SENOGRAPHE DMR MAMMOGRAPHY GANTRY SENOGRAPHE DMR CONTROL CONSOLE WITH RAD SHIELD SCREEN HYDRAULIC CHAIR	33 lbs 231 lbs 634 lbs 72 lbs	976 btu 4952 btu			

60		MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS 30 IN. W × 77 IN. H [762mm × 1956mm], CONTINGENT ON A 60 IN. [1524mm] CORRIDOR WIDTH					
61		MAIN DISCONNECT COUNTER TOP WITH BASE AND WALL CABINETS					

ENVIRONMENT

Ambient Operation Temperature: $50-104^{\circ}$ F $(10-40^{\circ}$ C) Allowable Temperature Change: 15° F $(10^{\circ}$ C) per hour

Humidity: 30-75 percent non condensing

Allowable Humidity Change: 10 percent per hour

Max - Altitude: 9500 ft. (29,000 M) above sea level